Coastal Plains Agri-Systems, LLC

Manuel d'instructions

CAS3

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- a. seconds of CO₂ gas flow, which will allow for an elevated CO₂ level of 75%.
- c. During the process, operator will notice sudden reflexes of the animal and tightening of the muscle which indicates shock. At this point, the operator should continuously administer CO₂ gas until the bird is completely still and appears to be lifeless and has stopped breathing.

The CAS-3 Stunning Device operates by replacing atmospheric oxygen with Carbon Dioxide. Carbon Dioxide is an odorless and tasteless gas that can cause physical harm if used improperly. Carefully read and follow all safety guidelines and warnings of this product below.

- d. Once the bird appears lifeless, remove the hood and perform the eye-touch test of the bird. If there is no reaction from the eye-touch, the bird will be considered as lifeless. Continue to monitor the bird for a few more seconds for reassurance of death.
 - e. In most cases, the time to death is 1 to 2 minutes, depending on the condition and age of the bird and the experience of the operator.
7. Cleaning and disinfection of your CAS-3:
- a. First, remove the seal by releasing the key drive clamp and removing the seal. Remove all organic material by hand washing before sanitizing.
 - b. Most of the parts of the CAS-3 may be power washed, or submerged in solutions except the pressure gauges and regulator.
 - c. Pressure regulators and gauges will need to be blown off with pressurized air and wiped clean by hand. Never submerge in a liquid solution.
 - d. Cleanup: First, normal ag detergents should be used for removal of organic material. Mild disinfectants such as quaternary ammonia products, iodine, and proper levels chlorine solution should be considered.
8. CO₂ first aid (CDC recommended first aid for CO₂ inhalation)
- a. If a person breathes large amounts of CO₂, move the exposed person to fresh air at once.
 - b. When breathing is difficult, properly trained personnel may assist the affected person by administering 100% oxygen.
 - c. If a person inhales large amounts of CO₂ and breathing has stopped, perform artificial respiration. Keep the affected person warm and at rest. Call 911 for medical assistance immediately.

Safety Precautions:

- Read and follow all label directions and operating instructions prior to the use of the CAS-3 euthanasia device.
- Do not operate the CAS-3 device in a non-ventilated structure. (OSHA 1910.146 (B).
- The CO₂ tank must be in off position when not in use.
- Prior to each use, audible and visual testing on all CO₂ tank fittings and hose fittings must be performed. This is a simple and quick test that can be achieved by pressurizing the system from the CO₂ tank on/off valve to the manual gas release valve on the hood.
 - a. Pressurize the system from the CO₂ tank's on/off valve to the manual gas release valve. Shut the tank valve off. Let stand for 1-2 minutes and observe the outlet pressure gauge. If the outlet pressure gauge pressure begins to drop, this indicates pressure loss.
 - b. The second step, mainly to conduct another type visual test would be to spray all fittings with a soap/water solution.
- After each cylinder exchange from empty to full, reset outlet pressure, as higher gas outlet pressure than 40 psi will result with a fully charged tank.
- Never place the CAS-3 hood over human orifices such as the nose or mouth.
- CO₂ can cause respiratory distress if misused. (See first aid in Section 8, Operating Instructions.)
- During heavy use and higher outlet pressure of 45 psi and above, operators should be aware of potential freezing of the regulators and CO₂ tanks. If frosting occurs, frostbite to the operator could result if touched. If freezing occurs, temporarily discontinue use of the CAS-3 until frosting has disappeared. Re-check the outlet pressure. Adjust outlet pressure if necessary to 40 psi or lower.
- Never store or transport CO₂ tanks inside vehicles or store in non-ventilated structures.
- CO₂ tanks must be properly secured to vehicle bed frames before attempting to transport.
- CO₂ tanks should be stored at temperatures below 125 degrees F / 51.7C.
- CO₂ first aid (CDC Recommended First Aid for CO₂ Inhalation)
 - a. If a person breathes large amounts of this chemical, moved the exposed person to fresh air at once. Other measures are usually unnecessary.
 - b. If breathing becomes difficult, move the exposed person to fresh air at once. When breathing is difficult, properly trained personnel may assist the affected person by administering 100 % oxygen. Keep the affected person warm and at rest.
 - c. If breathing has stopped, perform artificial respiration immediately. Keep the affected person warm and at rest. Get medical attention immediately by calling 911 for medical assistance.