

La qualité de la méthode pour euthanasier humainement un animal est une priorité pour l'utilisateur et pour le bien être de l'animal


OPÉRATION DU PISTOLET

Produit une contusion et trauma au cerveau, suffisant pour réduire une perte de conscience subite, puis la mort.

Le mouvement du piston et la cheville crée une énergie cinétique qui est transféré à l'animal par la cheville percutante



Zephyr



- Pistolet d'abattage à cheville percutante
- Au propane: 1 canne = 1000 coups
- 3 adaptateurs
> 40 lbs — 25-40 lbs — < 15 lbs

- **Dindes**
- **Poulets**
- **Oies**
- **Canards**
- **Porcelets**
- **Lapins**
- **< 20 lbs**

- insensibilité instantanée
- irréversible
- esthétique

➔ À utiliser avec compresseur à l'air




Euthanise humanly an animal today is the main priority for the users and for the Welfare of the animal.

HOW IT WORKS

Production of concussion and trauma to the cerebral hemisphere and brainstem, sufficient to induce sudden loss of consciousness and subsequent death (humane dispatch).

The kinetic energy is carried by the moving piston-bolt assembly and transferred to the subject by the bolt.



Zephyr



- Portable, captive Bolt working with propane can
- good for 1000 shots
- 3 adaptors
> 40 lbs — 25-40 lbs — < 15 lbs

- **Turkeys**
- **Chickens**
- **Geese**
- **Ducks**
- **Swine**
- **Rabbits**
- **< 20 lbs**

- instant insensibility
- irreversible
- esthetic

➔ To be use with a compressor 120 psi

